

EQUESTRIDE

EQUINE INJURY MANAGEMENT

SIZING GUIDE

The EqueStride™ tendon support device is available in three sizes - Small, Medium and Large.

To determine the correct size please provide the following two measurements:

1. The width of the fetlock joint at the widest point as viewed from the front
2. The width of the mid cannon bone, also as viewed from the front.
(See Example below)

In this example the width of the fetlock joint is 80mm (8CM) and the width of the mid cannon bone is 55mm (5.5CM)



1 FETLOCK JOINT WIDTH



2 CANNON BONE WIDTH

NOTE: Swelling of the limbs, (especially palmar swelling in the cannon area) can also affect the size requirement. Please let us know if there are any other abnormalities of the lower limb that could affect the sizing such as for example, a grossly enlarged tendon or a protruding splint bone. Where possible also send us photos of the limb to help us select the correct size device for the horse.